

Anti-Ikkε Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14164

Item Number: CM0010916

Introduction

Rabbit monoclonal antibody against IKKε (IKBKE), validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Recognizes N-terminal amino acids 1-100 of human IKKε (Q14164). Ideal for studying antiviral innate immunity and inflammatory signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	IKKε Rabbit mAb
Remarks/Alias	IKKE; IKK1; IKK-E; IKK-i; IKKε
Species	Human
GeneID (Human)	9641
GeneID	9641
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14164
Protein Weight	80kDa
Shipping	Ice bag