


Anti-Ikky Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q9Y6K9

Item Number: CM0010946



Introduction

High-affinity rabbit monoclonal antibody targeting human IKK γ (NF-kappa-B essential modulator). Validated for WB, IF/ICC, IHC-P, and ELISA. Detects endogenous levels in human, mouse, and rat. Essential for NF- κ B signaling pathway research.

[Learn More](#)

Parameter	Specification
Product Name	IKK γ Rabbit mAb
Remarks/Alias	IP, IP1, IP2, FIP3, IKKG, IPD2, NEMO, FIP-3, Fip3p, IMD33, SAIDX, AMCBX1, EDAID1, IKKAP1, ZC2HC9, IKK-gamma
Species	Human
GeneID (Human)	8517
GeneID	8517
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-419 of human IKK γ (NP_003630.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6K9
Protein Weight	48 kDa / 56 kDa / 37 kDa
Shipping	Ice bag