


Anti-Ikky Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y6K9

Item Number: CM0015157



Introduction

High-quality anti-IKKγ (NEMO) rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, and rat samples. Essential for NF-kappa-B signaling and IKK complex research.

[Learn More](#)

Parameter	Value
Product Name	IKKγ Rabbit pAb
Remarks/Alias	IP; IP1; IP2; FIP3; IKKG; IPD2; NEMO; FIP-3; Fip3p; IMD33; SAIDX; AMCBX1; EDAID1; IKKAP1; ZC2HC9; IKK-gamma; IKKγ
Species	Human
Gene ID	8517
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-390 of human IKKγ (NP_003630.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6K9
Protein Weight	48 kDa/56 kDa/37 kDa
Shipping	Ice bag