


Anti-Ikky Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Y6K9

Item Number: CM0013944



Introduction

Rabbit anti-IKKγ (NEMO) polyclonal antibody validated for Western Blot, Immunoprecipitation, and ELISA. Recognizes human and mouse NF-kappa-B essential modulator for studies in inflammatory signaling and disease research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] IKKγ Rabbit pAb
Remarks/Alias	IP; IP1; IP2; FIP3; IKKG; IPD2; NEMO; FIP-3; Fip3p; IMD33; SAIDX; AMCBX1; EDAID1; IKKAP1; ZC2HC9; IKK-gamma; Ky
Species	Human
GeneID (Human)	8517
GeneID	8517
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-419 of human IKKγ (NP_003630.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y6K9
Protein Weight	48 kDa/56 kDa/37 kDa
Shipping	Ice bag