


Anti-Ikk β Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O14920

Item Number: CM0007566



Introduction

High-quality rabbit monoclonal antibody targeting human IKK β (O14920) for Western blot, IHC-P, and ELISA applications. Ideal for NF- κ B signaling studies and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	IKK β Rabbit mAb
Remarks/Alias	IKK2; IKKB; IMD15; IMD15A; IMD15B; NFKB1KB; IKK-beta; IKK β
Species	Human
Gene ID (Human)	3551
Gene ID	3551
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human IKK β (O14920).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O14920
Protein Weight	29-87 kDa
Shipping	Ice bag