


Anti-Ikk β Rabbit Pab For Wb, If/Icc, Elisa - O14920

Item Number: CM0023125



Introduction

IKK β Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, and ELISA in human and mouse samples. It targets the IKK β kinase essential for NF-kB inflammatory signaling.

[Learn More](#)

Parameter	Specification
Product Name	IKK β Rabbit pAb
Remarks/Alias	IKK2, IKKB, IMD15, IMD15A, IMD15B, NFKB1KB, IKK-beta, IKK β
Species	Human
Gene ID (Human)	3551
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-669 of human IKK β (NP_001547.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14920
Protein Weight	29-87 kDa
Shipping	Ice bag