


Anti-Ikzf2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uks7

Item Number: CM0016932



Introduction

High-quality rabbit polyclonal antibody against human IKZF2 (Helios) validated for WB and ELISA. Cross-reacts with mouse and rat. Study Treg stability, IL-2 repression, and hearing development. Q9UKS7.

[Learn More](#)

Parameter	Value
Product Name	IKZF2 Rabbit pAb
Remarks/Alias	ANF1A2; HELIOS; ZNF1A2; ZNFN1A2; IKZF2
Species	Human
GeneID (Human)	22807
GeneID	22807
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-440 of human IKZF2 (NP_001072994.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKS7
Protein Weight	60kDa
Shipping	Ice bag