

Anti-Ikzf3 Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - Q9Ukt9

Item Number: CM0009986

Introduction

Anti-IKZF3 Monoclonal Antibody validated for WB, IP, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat IKZF3. Recombinant immunogen. Essential for lymphocyte differentiation, B-cell development, and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	IKZF3 Rabbit mAb
Remarks/Alias	AIO; IMD84; AIOLOS; ZNFN1A3; IKZF3
Species	Human
Gene ID (Human)	22806
Gene ID	22806
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-509 of human IKZF3 (Q9UKT9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKT9
Protein Weight	58 kDa
Shipping	Ice bag