

Anti-Ikzf3 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Ukt9

Item Number: CM0006863



Introduction

Rabbit monoclonal antibody targeting human IKZF3 (Aiolos) for WB, IF-P, IHC-P, and ELISA. Reacts with human and rat. Ideal for lymphocyte differentiation and B-cell biology research.

[Learn More](#)

Parameter	Value
Product Name	IKZF3 Rabbit mAb
Remarks/Alias	AIO; IMD84; AIOLOS; ZNFN1A3; IKZF3
Species	Human
Gene ID (Human)	22806
Gene ID	22806
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human IKZF3 (Q9UKT9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UKT9
Protein Weight	58kDa
Shipping	Ice bag