

Anti-Ikzf3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ukt9

Item Number: CM0026196

Introduction

Rabbit polyclonal antibody against IKZF3 (Aiolos), validated for WB, IF/ICC, and ELISA in human and mouse samples. Plays a key role in B-cell differentiation and apoptosis regulation.

[Learn More](#)

Parameter	Specification
Product Name	IKZF3 Rabbit pAb
Remarks/Alias	AIO; IMD84; AIOLOS; ZNFN1A3; IKZF3
Species	Human
GeneID (Human)	22806
GeneID	22806
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-450 of human IKZF3 (NP_036613.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UKT9
Protein Weight	58kDa
Shipping	Ice bag