


Anti-IL-1 Alpha Rabbit Polyclonal Antibody For Wb And Elisa - P01583

Item Number: CM0021663



Introduction

Rabbit polyclonal antibody against human IL-1 alpha, suitable for WB and ELISA. Validated for human, mouse, and rat samples. Immunogen corresponds to amino acids 112-271. Useful for studying inflammation and immune response.

[Learn More](#)

Parameter	Value
Product Name	IL-1 alpha Rabbit pAb
Remarks / Alias	IL1; IL-1A; IL1F1; IL1-ALPHA; IL-1 alpha
Species	Human
GeneID (Human)	3552
GeneID	3552
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 112-271 of human IL-1 alpha (NP_000566.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01583
Protein Weight	31kDa
Shipping	Ice bag