

Anti-IL-1 Beta Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P01584

Item Number: CM0015209



Introduction

Rabbit polyclonal antibody targeting human IL-1 beta (IL1 β , IL-1, IL1F2). Suitable for WB, IF/ICC, IF-P, and ELISA. Cross-reacts with human, mouse, and rat IL-1 β . For inflammation and cytokine research.

[Learn More](#)

Parameter	Value
Product Name	IL1 β Rabbit pAb
Remarks/Alias	IL-1; IL1F2; IL1beta; IL1-BETA; IL1 β
Species	Human
GeneID (Human)	3553
GeneID	3553
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-269 of human IL1 β (NP_000567.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01584
Protein Weight	31 kDa
Shipping	Ice bag