

Anti-IL-13 Rabbit Monoclonal Antibody For Fc, Elisa - P20109

Item Number: CM0006207

Introduction

Rabbit monoclonal antibody targeting mouse IL-13, validated for intracellular flow cytometry and ELISA. Specific to mouse IL-13, protein weight 14kDa. Suitable for allergy and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	IL-13 Rabbit mAb
Remarks/Alias	IL-13
Species	Mouse
GeneID (Human)	3596
GeneID	16163
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-131 of mouse IL-13 (NP_032381.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra), ELISA
Cross Reactivity	Mouse
SWISS	P20109
Protein Weight	14kDa
Shipping	Ice bag