


Anti-IL-13R α 1 (Cd213A1) Monoclonal Antibody For Flow Cytometry - P78552

Item Number: CM0008055



Introduction

-conjugated rabbit monoclonal antibody targeting human IL-13R α 1 (CD213a1). Validated for flow cytometry on human samples, supporting IL-13/IL-4 signaling studies in immunology and oncology.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD213a1/IL-13R α 1 mAb
Remarks/Alias	IL13RA1; CD213A1; CT19; IL-13Ra; NR4; interleukin-13 receptor subunit alpha-1
Species	Human
GeneID (Human)	3597
GeneID	3597
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 29-343 of human CD213a1/IL-13R α 1 (NP_001551.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P78552
Protein Weight	31kDa/48kDa
Shipping	Ice bag