

## Anti-IL-15Ra Rabbit Monoclonal Antibody For If/Icc, Fc, Elisa - Q60819

Item Number: CM0009581



### Introduction

High-affinity rabbit monoclonal antibody against mouse IL-15RA (CD215) for IF/ICC, FC, ELISA. Recombinant protein immunogen, MW 28kDa, cross-reacts with mouse. Suitable for studying IL-15 signaling, immune cell function.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | IL-15RA Rabbit mAb                                                                                                        |
| Remarks/Alias    | Il15ra                                                                                                                    |
| Species          | Mouse                                                                                                                     |
| GeneID (Human)   | 3601                                                                                                                      |
| GeneID           | 16169                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 33-205 of mouse CD215/IL-15Rα (NP_032384.1) |
| Source           | Rabbit                                                                                                                    |
| Category         | Monoclonal Antibodies                                                                                                     |
| Application      | IF/ICC, FC, ELISA                                                                                                         |
| Cross Reactivity | Mouse                                                                                                                     |
| SWISS            | Q60819                                                                                                                    |
| Protein Weight   | 28kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |