

Anti-IL-1Ra Rabbit Monoclonal Antibody For Wb And Elisa - P25085

Item Number: CM0006374

Introduction

Mouse interleukin-1 receptor antagonist (IL-1RA/IL-1ra) rabbit monoclonal antibody for WB and ELISA. Cross-reacts with human. Suitable for inflammation and cytokine research.

[Learn More](#)

Parameter	Value
Product Name	IL-1RA Rabbit mAb
Remarks/Alias	IL-1ra; F630041P17Rik; IL-1RA
Species	Mouse
GeneID (Human)	3557
GeneID	16181
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-178 of mouse IL-1RA (NP_001034790.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P25085
Protein Weight	20 kDa
Shipping	Ice bag