

Anti-IL-1Racp Monoclonal Antibody For Wb, Elisa - Q9Nph3

Item Number: CM0005967



Introduction

Monoclonal rabbit antibody targeting human IL-1 receptor accessory protein (IL-1RacP/IL1R3). Validated for WB, ELISA with cross-reactivity to mouse and rat. Ideal for IL-1, IL-36, and IL-33 signaling studies.

[Learn More](#)

Parameter	Value
Product Name	IL-1RacP Rabbit mAb
Remarks/Alias	IL1R3; C3orf13; IL-1RacP
Species	Human
GeneID (Human)	3556
GeneID	3556
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 391-470 of human IL-1RacP (NP_002173.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPH3
Protein Weight	65kDa
Shipping	Ice bag