


Anti-IL-1RacP Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nph3

Item Number: CM0030555



Introduction

Rabbit polyclonal antibody targeting human IL-1RacP (IL1R3), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for IL-1/IL-33/IL-36 signaling and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	IL-1RacP Rabbit pAb
Remarks/Alias	IL1R3; C3orf13; IL-1RacP
Species	Human
Gene ID (Human)	3556
Gene ID	3556
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-356 of human IL-1RacP (NP_001161402.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPH3
Protein Weight	65kDa
Shipping	Ice bag