

Anti-IL-21 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - Q9Hbe4

Item Number: CM0008642



Introduction

Anti-human IL-21 rabbit monoclonal antibody validated for WB, IF/ICC, IF-P, ELISA. Reacts with human and rat. Targets IL-21 cytokine involved in immunity, Tfh cell generation, and B-cell antibody production. Suitable for immunology research.

[Learn More](#)

Parameter	Value
Product Name	IL-21 Rabbit mAb
Remarks/Alias	Za11; IL-21; CVID11
Species	Human
Gene ID (Human)	59067
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-162 of human IL-21 (NP_068575.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9HBE4
Protein Weight	18kDa/19kDa
Shipping	Ice bag