

Anti-IL-27Ra/Wsx-1 Rabbit Monoclonal Antibody For Fc, Elisa - Q6Uwb1

Item Number: CM0005550



Introduction

Rabbit monoclonal antibody against human IL-27RA/WSX-1, validated for flow cytometry and ELISA with cross-reactivity to rat. Ideal for Th1 immune response and cytokine signaling research.

[Learn More](#)

Parameter	Value
Product Name	IL-27RA/WSX-1 Rabbit mAb
Remarks/Alias	CRL1; TCCR; WSX1; IL27R; IL-27RA; zcytor1; IL-27RA/WSX-1
Species	Human
GeneID (Human)	9466
GeneID	9466
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 33-283 of human IL-27RA/WSX-1, NP_004834.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q6UWB1
Protein Weight	69kDa
Shipping	Ice bag