


Anti-IL-33R/St2 Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P14719

Item Number: CM0004228



Introduction

Anti-IL-33R/ST2 rabbit monoclonal antibody validated for WB, IF/ICC, FC, ELISA. Recognizes mouse IL-33 receptor alpha chain (UniProt P14719) for innate and adaptive immunity research.

[Learn More](#)

Parameter	Value
Product Name	IL-33R/ST2 Rabbit mAb
Remarks/Alias	T1, St2, DER4, Ly84, ST2L, Fit-1, T1/ST2, St2-rs1
Species	Mouse
Gene ID (Human)	9173
Gene ID	17082
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-332 of mouse IL-33R/ST2 (NP_001020773.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Mouse
SWISS	P14719
Protein Weight	65kDa
Shipping	Ice bag