


Anti-IL-36 Gamma Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Nzh8

Item Number: CM0005060



Introduction

Recombinant rabbit monoclonal antibody targeting human IL-36 gamma (Q9NZH8). Validated for WB, IF/ICC, and ELISA applications. Ideal for studying skin inflammation, psoriasis, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	IL-36 Gamma Rabbit mAb
Remarks/Alias	IL36G; IL-1F9; IL-1H1; IL-1RP2; IL1E; IL1F9; IL1H1; IL1RP2; interleukin-36 gamma; IL-36 Gamma
Species	human
Gene ID (Human)	56300
Gene ID	56300
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-169 of human IL-36 Gamma (NP_062564.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9NZH8
Protein Weight	14kDa/18kDa
Shipping	Ice bag