

Anti-IL-6 Polyclonal Antibody For Wb And Elisa - P20607

Item Number: CM0025918

Introduction

IL-6 Rabbit pAb, validated for WB and ELISA applications. Reacts with human, mouse, and rat IL-6, an essential cytokine in immunity, inflammation, and metabolism.

[Learn More](#)

Parameter	Description
Product Name	IL-6 Rabbit pAb
Remarks/Alias	IL6; Ifnb2
Species	Rat
GeneID (Human)	3569
GeneID	24498
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-211 of rat IL6 (NP_036721.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20607
Protein Weight	24 kDa
Shipping	Ice bag