


Anti-IL-6 Rabbit Monoclonal Antibody For Wb And Elisa - P05231

Item Number: CM0003740



Introduction

High-specificity IL-6 rabbit monoclonal antibody validated for Western blot (WB) and ELISA, detecting human IL-6. Ideal for research on immune response, inflammation, and metabolic regulation. Recombinant immunogen ensures consistent performance.

[Learn More](#)

Parameter	Value
Product Name	IL-6 Rabbit mAb
Remarks/Alias	CDF; HGF; HSF; BSF2; IL-6; BSF-2; IFNB2; IFN-beta-2; IL6
Species	Human
Gene ID (Human)	3569
Gene ID	3569
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P05231
Protein Weight	21-24 kDa
Shipping	Ice bag