

Anti-IL-6 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Elisa - P20607

Item Number: CM0008214

Introduction

Rabbit monoclonal antibody against rat IL-6, validated for WB, IP, IF/ICC, ELISA. Cross-reacts with mouse IL-6. Suitable for studying acute phase response, inflammation, and immune regulation. Sourced via CamelBio.

[Learn More](#)

Parameter	Value
Product Name	IL-6 Rabbit mAb
Remarks/Alias	ILg6; Ifnb2
Species	Rat
GenelD (Human)	3569
GenelD	24498
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-211 of rat IL-6 (NP_036721.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P20607
Protein Weight	24 kDa
Shipping	Ice bag