


Anti-IL-6 Rabbit Polyclonal Antibody For Wb, Elisa - P05231

Item Number: CM0017479



Introduction

IL-6 Rabbit Polyclonal Antibody detects human, mouse, and rat IL-6 in WB and ELISA. Recognizes a key acute-phase cytokine involved in immune response, inflammation, and metabolism. Ideal for studying infection, tissue injury, and autoimmune diseases.

[Learn More](#)

Parameter	Value
Product Name	IL-6 Rabbit pAb
Remarks/Alias	CDF; HGF; HSF; BSF2; IL-6; BSF-2; IFNB2; IFN-beta-2; IL6
Species	Human
GeneID (Human)	3569
GeneID	3569
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-212 of human IL6 (NP_000591.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05231
Protein Weight	24kDa
Shipping	Ice bag