


Anti-IL-6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P05231

Item Number: CM0013608



Introduction

IL-6 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets cytokine IL-6 (P05231) involved in immunity, inflammation, and metabolism. Ideal for diagnostic and research assays.

[Learn More](#)

Parameter	Value
Product Name	IL-6 Rabbit pAb
Remarks/Alias	CDF; HGF; HSF; BSF2; IL-6; BSF-2; IFNB2; IFN-beta-2; IL6
Species	Human
Gene ID (Human)	3569
Gene ID	3569
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 30-212 of human IL6 (NP_000591.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05231
Protein Weight	24kDa
Shipping	Ice bag