

Anti-IL10 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P22301

Item Number: CM0024659



Introduction

Rabbit polyclonal antibody against IL10 (Human, Mouse, Rat) validated for WB, IHC-P, ELISA. Ideal for studying anti-inflammatory cytokines in immune regulation, inflammation, and disease research.

[Learn More](#)

Parameter	Value
Product Name	IL10 Rabbit pAb
Remarks/Alias	CSIF; TGIF; GVHDS; IL-10; IL10A; IL10
Species	Human
Gene ID (Human)	3586
Gene ID	3586
Immunogen	Synthetic peptide, A synthetic peptide corresponding to a sequence within amino acids 19-178 of human IL10 (NP_000563.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22301
Protein Weight	21kDa
Shipping	Ice bag