

Anti-IL10Rb Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q08334

Item Number: CM0017699



Introduction

High-quality rabbit polyclonal antibody targeting human IL10RB. Validated for WB, IF/ICC, and ELISA applications. Reacts with human, mouse, and rat. Ideal for cytokine signaling and antiviral defense research.

[Learn More](#)

Parameter	Value
Product Name	IL10RB Rabbit pAb
Remarks/Alias	CRFB4; IBD25; CRF2-4; D21S58; D21S66; CDW210B; IL-10R2; IL-10RB; IL10RB
Species	Human
GeneID (human)	3588
GeneID	3588
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-220 of human IL10RB (NP_000619.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08334
Protein Weight	37kDa
Shipping	Ice bag