


Anti-IL11Ra Rabbit Mab For Wb - Q14626

Item Number: CM0011569



Introduction

Recombinant rabbit monoclonal antibody targeting human IL11RA (Q14626). Validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for cancer and bone development studies.

[Learn More](#)

Parameter	Value
Product Name	IL11RA Rabbit mAb
Remarks/Alias	CRSDA; IL11RA
Species	Human
Gene ID (Human)	3590
Gene ID	3590
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14626
Protein Weight	45kDa
Shipping	Ice bag