


## Anti-IL12 Rabbit Polyclonal Antibody For Wb, Elisa - P29459 / P29460

Item Number: CM0027314



### Introduction

Rabbit polyclonal antibody against human IL12 (p35/IL12A and p40/IL12B), suitable for WB and ELISA. Useful for studying IL-12 cytokine signaling, Th1 differentiation, and immune response mechanisms.

[Learn More](#)

| Parameter        | Value                                                                                                                                                  |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | IL12 Rabbit pAb                                                                                                                                        |
| Remarks/Alias    | P35; CLMF; NFSK; NKSF1; IL-12A                                                                                                                         |
| Species          | Human                                                                                                                                                  |
| GeneID (Human)   | 3592/3593                                                                                                                                              |
| GeneID           | 3592/3593                                                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 23-328 of human IL12B/23-219 of human IL12A (NP_002178.2/NP_001384921.1) |
| Source           | Rabbit                                                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                                                  |
| Application      | WB, ELISA                                                                                                                                              |
| Cross Reactivity | Human                                                                                                                                                  |
| SWISS            | P29459/P29460                                                                                                                                          |
| Protein Weight   | 25kDa/37kDa                                                                                                                                            |
| Shipping         | Ice bag                                                                                                                                                |