


Anti-IL12A Polyclonal Antibody For Wb, Elisa - P29459

Item Number: CM0024075



Introduction

IL12A Rabbit pAb (P29459) is a polyclonal antibody raised against a synthetic peptide within amino acids 50-150 of human IL12A. Suitable for WB and ELISA with cross-reactivity to human, mouse, and rat IL12A.

[Learn More](#)

Parameter	Value
Product Name	IL12A Rabbit pAb
Remarks/Alias	P35; CLMF; NFSK; NKSF1; IL-12A; IL12A
Species	Human
GeneID (Human)	3592
GeneID	3592
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human IL12A (NP_000873.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29459
Protein Weight	25kDa
Shipping	Ice bag