


Anti-IL12B Rabbit Monoclonal Antibody For Wb, Ihc-P, Fc, Elisa - P43432

Item Number: CM0004438



Introduction

Rabbit monoclonal antibody targeting mouse IL12B (p40), validated for WB, IHC-P, FC (intra), and ELISA. Cross-reacts with human. Ideal for IL-12/IL-23 signaling research.

[Learn More](#)

Parameter	Value
Product Name	IL12B Rabbit mAb
Remarks/Alias	p40; IL-12b; IL12p40; IL-12p40; IL12B
Species	Mouse
Gene ID (Human)	3593
Gene ID	16160
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-335 of mouse IL12B(NP_001290173.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC (intra), ELISA
Cross Reactivity	Human, Mouse
SWISS	P43432
Protein Weight	38kDa
Shipping	Ice bag