

Anti-IL12B Rabbit Polyclonal Antibody For Wb, Elisa - P29460

Item Number: CM0014922



Introduction

IL12B Rabbit pAb validated for WB and ELISA. Recognizes human and mouse IL-12 subunit beta, a cytokine involved in T and NK cell activation and IFN-gamma production. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	IL12B Rabbit pAb
Remarks/Alias	CLMF; NKSF; CLMF2; IMD28; IMD29; NKSF2; IL-12B; IL12B
Species	Human
Gene ID (Human)	3593
Gene ID	3593
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-328 of human IL12B (NP_002178.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P29460
Protein Weight	37kDa
Shipping	Ice bag