


## Anti-IL12B Rabbit Polyclonal Antibody For Wb, Elisa - P43432

Item Number: CM0020798



### Introduction

Mouse IL12B rabbit polyclonal antibody for WB and ELISA, cross-reacting with human. Targets the p40 subunit of IL-12 and IL-23, key cytokines in T/NK cell activation and IFN-gamma production. Research use only.

[Learn More](#)

| Parameter        | Value                                                                                                                                         |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | IL12B Rabbit pAb                                                                                                                              |
| Remarks/Alias    | p40; IL-12b; IL12p40; IL-12p40; IL12B                                                                                                         |
| Species          | Mouse                                                                                                                                         |
| GeneID (Human)   | 3593                                                                                                                                          |
| GeneID           | 16160                                                                                                                                         |
| Immunogen        | Recombinant protein<br>Recombinant fusion protein containing a sequence corresponding to amino acids 23-335 of mouse IL12B. (NP_001290173.1). |
| Source           | Rabbit                                                                                                                                        |
| Category         | Polyclonal Antibodies                                                                                                                         |
| Application      | WB, ELISA                                                                                                                                     |
| Cross Reactivity | Human                                                                                                                                         |
| SWISS            | P43432                                                                                                                                        |
| Protein Weight   | 38kDa                                                                                                                                         |
| Shipping         | Ice bag                                                                                                                                       |