


Anti-IL12Rb1 Rabbit Polyclonal Antibody For Wb, Elisa - P42701

Item Number: CM0023627



Introduction

High-quality rabbit polyclonal antibody targeting human IL12RB1, suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Useful for studying IL-12/IL-23 signaling and immune response mechanisms.

[Learn More](#)

Parameter	Value
Product Name	IL12RB1 Rabbit pAb
Remarks/Alias	CD212; IMD30; IL12RB; IL-12R-BETA1; IL12RB1
Species	Human
Gene ID (Human)	3594
Gene ID	3594
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-305 of human IL12RB1 (NP_714912.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42701
Protein Weight	73kDa
Shipping	Ice bag