


Anti-IL12Rb2 Polyclonal Antibody For Wb, Elisa - Q99665

Item Number: CM0013069



Introduction

Human IL12RB2 rabbit polyclonal antibody for WB and ELISA, cross-reactive with mouse and rat. Ideal for studying IL-12 receptor signaling in T-cell and NK cell biology.

[Learn More](#)

Parameter	Value
Product Name	IL12RB2 Rabbit pAb
Remarks/Alias	IL12RB2
Species	Human
Gene ID (human)	3595
Gene ID	3595
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 420-620 of human IL12RB2 (NP_001550.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q99665
Protein Weight	97kDa
Shipping	Ice bag