


Anti-IL13Ra1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P78552

Item Number: CM0009386



Introduction

High-quality rabbit monoclonal antibody targeting human IL13RA1 (CD213A1, NR4, CT19). Validated for Western blot, IHC, and ELISA. Explore IL-13/IL-4 signaling and cancer/testis antigen biology.

[Learn More](#)

Parameter	Value
Product Name	IL13RA1 Rabbit mAb
Remarks/Alias	NR4; CT19; CD213A1; IL-13Ra; IL13RA1
Species	Human
GeneID (Human)	3597
GeneID	3597
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human IL13RA1 (P78552).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P78552
Protein Weight	49kDa
Shipping	Ice bag