


Anti-IL13Ra1 Polyclonal Antibody For Wb, Elisa - P78552

Item Number: CM0027921



Introduction

Rabbit polyclonal antibody targeting interleukin-13 receptor subunit alpha-1 (IL13RA1, CD213A1). Suitable for WB and ELISA in human and mouse. Recombinant immunogen detects ~49kDa protein. Ideal for IL13/IL4 signaling and inflammatory disease research.

[Learn More](#)

Parameter	Value
Product Name	IL13RA1 Rabbit pAb
Remarks/Alias	NR4; CT19; CD213A1; IL-13Ra; IL13RA1
Species	Human
GeneID (human)	3597
GeneID	3597
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 29-343 of human IL13RA1 (NP_001551.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P78552
Protein Weight	49kDa
Shipping	Ice bag