

Anti-IL13Ra2 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q14627

Item Number: CM0021274



Introduction

Anti-IL13RA2 Rabbit Polyclonal Antibody validated for WB, IF/ICC, IHC-P, ELISA. Detects human, mouse, rat IL13RA2 protein. Relevant for IL-13 signaling and immune modulation studies.

[Learn More](#)

Parameter	Value
Product Name	IL13RA2 Rabbit pAb
Remarks/Alias	CT19; IL-13R; IL13BP; CD213A2; IL13RA2
Species	Human
Gene ID	3598 (Human)
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14627
Protein Weight	44kDa
Shipping	Ice bag