

Anti-IL13Ra2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14627

Item Number: CM0019592



Introduction

Rabbit polyclonal antibody targeting human IL13RA2, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying IL-13 signaling, immune regulation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	IL13RA2 Rabbit pAb
Remarks/Alias	CT19; IL-13R; IL13BP; CD213A2; IL13RA2
Species	Human
Gene ID (Human)	3598
Gene ID	3598
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-343 of human IL13RA2 (NP_000631.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14627
Protein Weight	44kDa
Shipping	Ice bag