

Anti-IL15Ra Rabbit Polyclonal Antibody For Wb, Elisa - Q60819

Item Number: CM0020499



Introduction

Rabbit polyclonal antibody targeting mouse interleukin-15 receptor subunit alpha (IL-15RA). Validated for WB and ELISA applications. Suitable for studies on IL-15 signaling, immune regulation, and cytokine receptor biology.

[Learn More](#)

Parameter	Value
Product Name	IL15ra Rabbit pAb
Remarks/Alias	IL-15RA; IL15ra
Species	Mouse
GeneID (Human)	3601
GeneID	16169
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 33-105 of mouse IL15ra (NP_032384.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q60819
Protein Weight	28kDa
Shipping	Ice bag