

Anti-IL17A Rabbit Monoclonal Antibody For Wb, Elisa - Q62386

Item Number: CM0007333

Introduction

Rabbit monoclonal antibody targeting mouse IL-17A, validated for WB and ELISA. Recognizes a 17 kDa protein involved in Th17 immunity and inflammation. Ideal for immune response and host defense research.

[Learn More](#)

Parameter	Value
Product Name	IL17A Rabbit mAb
Remarks/Alias	Il17; CtlA8; IL-17; CtlA-8; IL-17A; IL17A
Species	Mouse
Gene ID	16171
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-158 of mouse IL17A. (NP_034682.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q62386
Protein Weight	17kDa
Shipping	Ice bag