

Anti-IL17C Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9P0M4

Item Number: CM0014952

Introduction

This IL17C rabbit polyclonal antibody, validated for WB, IF/ICC, and ELISA, recognizes human, mouse, and rat IL-17C. Ideal for studying epithelial innate immunity and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	IL17C Rabbit pAb
Remarks/Alias	CX2; IL-17C; IL17C
Species	Human
Gene ID (Human)	27189
Gene ID	27189
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-197 of human IL17C (NP_037410.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P0M4
Protein Weight	22kDa
Shipping	Ice bag