

Anti-IL17F Rabbit Polyclonal Antibody For Wb, Elisa - Q96Pd4

Item Number: CM0030133



Introduction

Rabbit polyclonal antibody against human IL17F, suitable for WB and ELISA. Reacts with human and mouse; targets a key Th17 cytokine in innate and adaptive immunity.

[Learn More](#)

Parameter	Value
Product Name	IL17F Rabbit pAb
Remarks/Alias	ML1; ML-1; IL17A; CANDF6; IL-17F
Species	Human
GeneID (Human)	112744
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-163 of human IL17F (NP_443104.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96PD4
Protein Weight	18kDa
Shipping	Ice bag