

Anti-IL17Ra Rabbit Polyclonal Antibody For Wb, Elisa - Q96F46

Item Number: CM0015278



Introduction

Rabbit polyclonal antibody against human IL17RA (IL-17 receptor A, CD217), validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for studying IL-17 signaling, inflammation, and host defense.

[Learn More](#)

Parameter	Value
Product Name	IL17RA Rabbit pAb
Remarks/Alias	CD217; IL17R; IMD51; CANDF5; CDw217; IL-17RA; hIL-17R; IL17RA
Species	Human
GeneID (human)	23765
GeneID	23765
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-320 of human IL17RA (NP_055154.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96F46
Protein Weight	96 kDa
Shipping	Ice bag