


Anti-IL17Rb Rabbit Polyclonal Antibody For Wb & Elisa - Q9Nrm6

Item Number: CM0017854



Introduction

Rabbit polyclonal antibody specific for human IL17RB, validated for WB and ELISA. Exhibits cross-reactivity with mouse and rat. Target is a receptor for IL17B and IL17E, involved in hematopoietic cell growth/differentiation.

[Learn More](#)

Parameter	Value
Product Name	IL17RB Rabbit pAb
Remarks/Alias	CRL4; EVI27; IL17BR; IL17RH1; IL17RB
Species	Human
Gene ID (Human)	55540
Gene ID	55540
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-292 of human IL17RB (NP_061195.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRM6
Protein Weight	56kDa
Shipping	Ice bag