

Anti-IL17Rd Polyclonal Antibody For Wb And Elisa - Q8Nfm7

Item Number: CM0017827

Introduction

Rabbit polyclonal antibody targeting human IL17RD (SEF, IL17RLM), validated for WB and ELISA. Detects human, mouse, and rat IL17RD. Suitable for studying FGF signaling, Ras-MAPK pathway inhibition, and receptor localization.

[Learn More](#)

Parameter	Value
Product Name	IL17RD Rabbit pAb
Remarks/Alias	SEF; HH18; IL-17RD; IL17RLM; IL17RD
Species	Human
GeneID (Human)	54756
GeneID	54756
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-300 of human IL17RD (NP_060033.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NFM7
Protein Weight	82kDa
Shipping	Ice bag