

Anti-IL17Re Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nfr9

Item Number: CM0016216

Introduction

Rabbit polyclonal antibody against human IL17RE, validated for WB and ELISA applications. Recognizes human and mouse IL17RE (UniProt Q8NFR9), a receptor for IL17C involved in innate immunity.

[Learn More](#)

Parameter	Value
Product Name	IL17RE Rabbit pAb
Remarks/Alias	IL17RE
Species	Human
GeneID	132014
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-260 of human IL17RE (NP_705613.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NFR9
Protein Weight	75kDa
Shipping	Ice bag