

Anti-IL18 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P70380

Item Number: CM0021736

Introduction

Rabbit polyclonal antibody against mouse IL18, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Targets a pro-inflammatory cytokine involved in Th1 and NK cell responses.

[Learn More](#)

Parameter	Value
Product Name	IL18 Rabbit pAb
Remarks/Alias	Ig1f; IL-18; IL18
Species	Mouse
Gene ID (Human)	3606
Gene ID	16173
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P70380
Protein Weight	22kDa
Shipping	Ice bag